

C3orf57-chromosome 3 open reading frame 57 Gene

Reference: PTXBC065208

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf57; ADMP; SSSPTB; serine palmitoyltransferase small subunit B; androgen down regulated in mouse prostate; likely ortholog of androgen down regulated gene expressed in mouse prostate; small subunit of serine palmitoyltransferase B		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 57
GENE ID	165679
ORIGIN SPECIES	Human
ACCESSIONS	BC065208
VECTOR	pDONR223
NCBI SYMBOL	C3orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattgagcgctgtgaaggaatattctcctgctctactatcaataaccaatcattagctgctgtgctgttttagaccctggagcgcctatgtttaaccacatcttaaccattatgctatggtgtatatactgctatgctttatccaatccacattgcctggcttgggaattttctc
aaaaatatgtgatatcacagtacaatttctaattga