



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF323-zinc finger protein 323 Gene

Reference: PTXBC065241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF323; ZNF20-Lp; ZNF310P; zinc finger and SCAN domain-containing protein 31; zinc finger protein 323; zinc finger and SCAN domain containing 31		

For research use only.

GENE DESCRIPTION	zinc finger protein 323
GENE ID	64288
ORIGIN SPECIES	Human
ACCESSIONS	BC065241
VECTOR	pDONR223
NCBI SYMBOL	ZNF323
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcacacacgcacgatgatctcttttccctccaccaattcaagaacaagatgtgaaagtatactcgtgagaacacaggagttggcactcaaaagcaagaacctttaaagaagaatggaacattgggggatagcaaacctcaaaagagatgtatcttggattctaagtacagagaacctgttaa
acgagacacgaaggcgaaaaaacgacgagcgcacattccactggagagagacgccacaggctcaatgaatgtgggaaaagctcactaagatgtcagctactcattgagcaccagagaaatccacactgggggaagcccatatgaatgtgaaagatgtgggaagcgcttcagccggaggtc
aagctcgaatgaacatcgccggagccacactggagagaacctatcaatgtaagagatgtgggaaaagccttcagtgccagcaatggcctcactcagacagaagaatcaacacagaggggaaaaaacatagatgaaagtgtgtgggaaagccttccctcagctcatgcctgttca
gcatcagagagatcacagagagaagcgcctactcgtgtcgtgtgtggcaaaagccttcattcagatgcagagcgtttccagcatcctcaggctccacactggtgagaacctcatcgtgcagtcagtcgagtaaccttttagtaagcggacactcttaagaacatcagaaaaacca
cactgacagagagacataaa

More info: Please contact us for any question: contact@proteogenix.fr