



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GTF2H3-general transcription factor IIH, polypeptide 3, 34kDa Gene

Reference: PTXBC065250

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTF2; P34; TFB4; TFIIH; general transcription factor IIH subunit 3; TFIIH basal transcription factor complex p34 subunit; basic transcription factor 2 34 kDa subunit; general transcription factor IIH, polypeptide 3, 34kDa		

For research use only.

GENE DESCRIPTION	general transcription factor IIH, polypeptide 3, 34kDa
GENE ID	2967
ORIGIN SPECIES	Human
ACCESSIONS	BC065250
VECTOR	pDONR223
NCBI SYMBOL	GTF2H3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttttcagacagaagatgaattgaattcttctgttatgttagttgtatgcccaaccaatttgggtgggaagacgaacattaaaggaatctcagttacttattccaatgcatagatccgtgtatggtgtctgggaaatcgcattttcatgaatcgttccacaacattcgtgtatagcaagtcac
attcaagaagaacgcgattctatattcctgggaagaatggcagacttggagacttctcggagaccctggcaaccctcctgaatttaacctctgggagtaaaagtggaataatcgaacttttaacctcagcaaatgaagtattgttgaagaagattaaagatctaagaccaaaagtgcataaaag
gttcacatcacacagatctctctgcgcaggatccctgcgccaagcccttctctacatcattagaatgaacacagggaagttaaagacacatcaggaaatgaataacaggatattggtgtatgaagctgcagaagacagctgtgttcgagfatagcaatgaatgcatttccagcagacaagaac
agaatttttattgattgagcctgtgtttttagactccctgagcctccccaacagcggggagctgacctgaagctgacctcagatcctctctctgcgatttctgtggtgtgttctcccgatcaagatcagagatcctgatattccctcccccgggttcattgtga
ctacagcgtctgtctctctctcagaactcattgaaattgtatgtctctctctgttctgttgcattatctcaatttccagcccatttctacatctctgcagacagcccttaaaatttctctccctcagctctgaagccagaagaagaactgaattgttcttccctga

More info: Please contact us for any question: contact@proteogenix.fr