

MRPL54-mitochondrial ribosomal protein L54 Gene

Reference: PTXBC065273

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L54mt; MRP-L54; 39S ribosomal protein L54, mitochondrial; mitochondrial ribosomal protein L54		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L54
GENE ID	116541
ORIGIN SPECIES	Human
ACCESSIONS	BC065273
VECTOR	pDONR223
NCBI SYMBOL	MRPL54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaccaaaccgttttggggctacccggacgtgggcccggctggggggcctggagctcctaaacccgccacttcggaagactcctggcccggattatgccaagaaccagttatgaagggggccaaatcgggaaaaggtgcagtgaccagcgaggccctcaaggacc
ccgacgtatgcacagatcctgtccagctcaccacatgcccattggcgtcaacatctacaaggaaggcaggatgtacccctgaaaccggatgctgagtagcctgaatggctgtttgagatgaacttgggtcccccagaccctggaggagctggaccccgagagccgggagtagtgg
cggcggctgcggaaacagaacatctggcgccacaaccggctgagcaagaacaagaggtgtag