

HSPA14-heat shock 70kDa protein 14 Gene

Reference: PTXBC065281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSP70-4; HSP70L1; heat shock 70 kDa protein 14; HSP70-like protein 1; heat shock 70kDa protein 14 isoform 1 variant 3; heat shock protein HSP60; heat shock protein hsp70-related protein; heat shock protein family A (Hsp70) member 14		

For research use only.

GENE DESCRIPTION	heat shock 70kDa protein 14
GENE ID	51182
ORIGIN SPECIES	Human
ACCESSIONS	BC065281
VECTOR	pDONR223
NCBI SYMBOL	HSPA14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccatcgaggttcacgtggcgtgcacctcagcctgtgtgcccgtctataaggatggcgggctgtgtgttgcaaatgatgccggtgaccgagttactcagctgtgtgtcttactcagaaatgaagagattgttgattggcagcaaaacaaagtagaataagaaatttcaa
atacagtaataagtaagcagatcctgggcagaagccagaatgcgtcttgacctggccttctcagcaattatccctga