

BCAP31-B-cell receptor-associated protein 31 Gene

Reference: PTXBC065292

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6C6-AG; BAP31; CDM; DDCH; DXS1357E; B-cell receptor-associated protein 31; 6C6-AG tumor-associated antigen; BCR-associated protein Bap31; p28 Bap31		

For research use only.

GENE DESCRIPTION	B-cell receptor-associated protein 31
GENE ID	10134
ORIGIN SPECIES	Human
ACCESSIONS	BC065292
VECTOR	pDONR223
NCBI SYMBOL	BCAP31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcgcagtggaactgcagttgccacctctctatgcggagggtctttgtgtgtgtctctctcattcccttcattctcctaaagatggcagaagatttcaagtcctccgctgggtgagttgttagtgcctatggcaacacctctttgtgtgtctcattgcctctgtgtgtgtgtcatcga
tgccgtgctgcgaaatcggaagtatgatgatgtgacggaaaaggtgaacctccagaacaatccggggccatggagcactccacatgaagctttccgtgccagaggaaatctctacattgctgcttttctctgtctgtccttctctgttagacgcctggtgactctcatttcgacgcagg
ccacgctgctggcctccaatgaagcctttaaaaagcaggcggagagtgctagtgaggcggccaagaagtacatggaggagaatgaccagctcaagaaggagctgctgttgacggaggaagttggatgtcgggaatgctgaggtgaagttggaggaagagaacaggagcctgaag
gctgacctgcagaagctaaaggacgagctggccagcactaagcaaaactagagaagctgaaaaccaggttctggccatgcgggaagcagctgagggcctaccaaggagtagcaccgctgtgctgaggagcacgcaagctgaggctgcagtagatggtccatggacaaga
aggaaagagtaa