

CCDC25-coiled-coil domain containing 25 Gene

Reference: PTXBC065295

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 25; coiled-coil domain containing 25		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 25
GENE ID	55246
ORIGIN SPECIES	Human
ACCESSIONS	BC065295
VECTOR	pDONR223
NCBI SYMBOL	CCDC25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgttctactcaccagcagcagcgttaattcatctgcctacactatttacatgggaaagataaatgaagatctgatcaagcatggctggcctgaagatatctgtttcatgtggacaaactctctcggtcatgtataccttcgattacataaggagagaatagaaagacat
cccaagggaagtgtgactgtgcccacctgtggaaggccaatgacattcaaggctgcaagatgaacaacgttaattgtgtatatacgcctgtgctaacctgaagaaaacagctgacatggatgtgggcagataggcttcacagcagaaggatgtaaaattgtgacagtggag
aagaaagtaaatgagatcctgaaccgattgaaaagacaaagtcgagcgggtccagacctagcagcagagaagaatgcagagatctgaagagagggaatgagaaaaagcccaattcagggaatgaaaagagagaaaaaagagaatgaagaagagggaatggatg
aacttaggagctattcatcactaatgaaagtggaaaatgtcttcaatcaggatggcaatgattcagatgaattcatgtaa