



Reference: PTXBC065300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIH003; UBP; ubiquitin carboxyl-terminal hydrolase 3; deubiquitinating enzyme 3; ubiquitin thioesterase 3; ubiquitin thiolesterase 3; ubiquitin-specific-processing protease 3; ubiquitin specific peptidase 3		

For research use only.

GENE DESCRIPTION	ubiquitin specific peptidase 3
GENE ID	9960
ORIGIN SPECIES	Human
ACCESSIONS	BC065300
VECTOR	pDONR223
NCBI SYMBOL	USP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]