

LOC100129534-small nuclear ribonucleoprotein polypeptide N pseudogene Gene

Reference: PTXBC065369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	small nuclear ribonucleoprotein polypeptide N pseudogene
GENE ID	100129534
ORIGIN SPECIES	Human
ACCESSIONS	BC065369
VECTOR	pDONR223
NCBI SYMBOL	LOC100129534
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctactggtgtccctcgagcagagggaagcccaattggtggcccggtgggtggcctaccaggttcagaagctccaagtttaactgtggacattggaatcggggaacggcagtcgtgacctgggcaagtagcaagatgctgccgcacacggaccacagaatgagacgtat
cctacaagatggccgagtgcttcattggcaccctcaaggcttttgacaagcatatgaatgatccctgtgattgtgacgagttcagaaaatgaaactaaagaatcggaagcaaccagcacgtgaaagaaagcgggtttcggggccctccgcttctaaagaagcttggggggccaca
ctcactgccattgctcagagcttggttacacggccggccacacagggccacacagagacaagaagcgcaatgtgctggcaaaactcttaccggaggagaaggggagaatggctgtcttggccaatagcccggtcacgaacgtaccctgctggggcccccgaagtgcaccaaccg
actcgaaggtccagagagggtgctggggcaaggacgtgagttccagccaccctgcaggaggagcccaggaccctgcacccctga