

CXorf40A-chromosome X open reading frame 40A Gene

Reference: PTXBC065539

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CXorf40A; CXorf40; EOLA1; FLJ16423; endothelial-overexpressed lipopolysaccharide-associated factor 1; chromosome X open reading frame 40A		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 40A
GENE ID	91966
ORIGIN SPECIES	Human
ACCESSIONS	BC065539
VECTOR	pDONR223
NCBI SYMBOL	CXorf40A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtttgctgcctctcctccggcagccttatgctgcttcttaaatggaatcaagactgtggagacgcgtggcgtccctgtgagcagccagcggaaactgtaccatgccgtccacattgtcacagggaactgggaaggcgatgcctggcgggagctgtgtggagagac
tcgggatgactctgtcagattcaggcctgtcaggaaggggaaaagtgtgctgaggagtgatagcgggactcgttgacattgggaaactttgcaatgccccgaagacttaactcccgaatgaggtgttggaactagaaaatcaagctgtactgaccaacctgaagcagaagtacgtg
actgtgatttcaaacccagggtggttactggagccatcctaggaaggaggcaaggatgtattccaggtagacatccagagcacctgatcccttggggcatgaagtgtga