



Reference: PTXBC065540

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APP-1; APP1; PABP4; iPABP; polyadenylate-binding protein 4; PABP-4; activated-platelet protein 1; inducible poly(A)-binding protein; poly(A) binding protein, cytoplasmic 4 (inducible form); poly(A)-binding protein 4; poly(A) binding protein cytoplasmic 4		

GENE DESCRIPTION	poly(A) binding protein, cytoplasmic 4 (inducible form)
GENE ID	8761
ORIGIN SPECIES	Human
ACCESSIONS	BC065540
VECTOR	pDONR223
NCBI SYMBOL	PABPC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]