

RD3-retinal degeneration 3 Gene

Reference: PTXBC065541

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein RD3; C1orf36; LCA12; retinal degeneration protein 3; retinal degeneration 3		
For research use only.			

GENE DESCRIPTION	retinal degeneration 3
GENE ID	343035
ORIGIN SPECIES	Human
ACCESSIONS	BC065541
VECTOR	pDONR223
NCBI SYMBOL	RD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctcatctatgcttcggtggaacgaggcccatcccgctgtccaccaggagccctgctgagatggtgctggagacgcttatgatggagctgacggggcagatgcgagaggctgagaggcagcagcggggagcgcagcaatgcggtcagaaaggctgcaccggtgtggac
tacagctggctggccagcacaccccggtcccatatgacctcagcccaattgagcgggtgacagctggaagatgtctgctgtaagatccaccatcctattgtggcctgctatcctcagggtccggcagctgctggcggagcaggagcccgagggtcaggagggtgtccagctcttccgct
cggctgctgcaggaggtcctggagaggatgaagcaggagaggagcccaagctgacgcccagctggagcctgctggcccccggcagcctggccacctcaagaccgcgcgcatctcgccttcgccagcgacatcaggacctctccgaggacgtggagcgggacaca
ccgccccactgcggtcctggagcatgcccgaattccggcgcccaagccgactga