



Reference: PTXBC065555

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TOMM70A; Tom70; mitochondrial import receptor subunit TOM70; mitochondrial precursor proteins import receptor; translocase of outer membrane 70 kDa subunit; translocase of outer mitochondrial membrane 70 homolog A; translocase of outer mitochondrial membr		

For research use only.

GENE DESCRIPTION	translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae)
GENE ID	9868
ORIGIN SPECIES	Human
ACCESSIONS	BC065555
VECTOR	pDONR223
NCBI SYMBOL	TOMM70A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr