

MGC70863-similar to RPL23AP7 protein Gene

Reference: PTXBC065556

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	similar to RPL23AP7 protein
GENE ID	284942
ORIGIN SPECIES	Human
ACCESSIONS	BC065556
VECTOR	pDONR223
NCBI SYMBOL	MGC70863
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcactcacttcaggcggccaagactcgcactccgaggcagcccagatcctcggagagcaccgccagagaacaagcttggccactatgctatcatcaagttccgctggccactgagtcggccgtgaagaagatagaagaaacaacacgcttgttctcactgtggat
gttaaagccaacaagcaccagatcagacaggctgtgaagaagctctatgacagtgatgtggccaaggtcaccacctgattgtctgataaggagaacaagcatatgttcgacttgcctctgattatgatgctttcgtgttaacaaaattgggatcacctaaactgagtcacagctgcta
a