

TP53I11-tumor protein p53 inducible protein 11 Gene

Reference: PTXBC065557

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PIG11; tumor protein p53-inducible protein 11; p53-induced gene 11 protein; tumor protein p53 inducible protein 11		

For research use only.

GENE DESCRIPTION	tumor protein p53 inducible protein 11
GENE ID	9537
ORIGIN SPECIES	Human
ACCESSIONS	BC065557
VECTOR	pDONR223
NCBI SYMBOL	TP53I11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacacccggaagtggcgagcccaagccctggggacaggtgtagggagaaaagcagcccgagcctcagactcgctccatcctggcatagagtgggggatggctggagggtgtctataggtacagcccgtctggctgtgccaggtggccctgccagggtctcacc
cctgtccaccctgtgctggtgtgctgctgacccagatagcaacatggcctgtaccagcagagtggtggcaccacctggttacagcggatgccccgagactctgctgtgtaaacgtggcagagcagaatgggaggtggaccctgaggaggccctctctctgcatctgtctct
tgctacctaa