



Reference: PTXBC065566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SVAS; WBS; tropoelastin		
For research use only.			

GENE DESCRIPTION	elastin
GENE ID	2006
ORIGIN SPECIES	Human
ACCESSIONS	BC065566
VECTOR	pDONR223
NCBI SYMBOL	ELN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1