

## ASB5-ankyrin repeat and SOCS box-containing 5 Gene

Reference: PTXBC065710

|                |                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ASB-5; ankyrin repeat and SOCS box protein 5; SOCS box protein ASB-5; ankyrin repeat and SOCS box containing 5 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ankyrin repeat and SOCS box-containing 5                                                          |
| <b>GENE ID</b>          | 140458                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC065710                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ASB5                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtcgggttagaagaaatcggccgtttgctcaacaattatccaatgtctactttacaatacttgcgtgttctgttttaagctttttgtgaaatcagccttgccatcctcagtcatttctacatagtgaaggcaaccgcaaggagcggaaggatagcagctgaattttatgagtaaccaag  
gacaagggttcctggcgagatcgcaccactacatgaagcagcaagtcaggtcgccttctgtctgagaacattattatcacagggttataatgtaaatgcagtaaccttagaccatgtcaccccatgtcacgaagcctgtgagatcacgtggcatgtgccagaactctgtggaagc  
aggagctaatgtaaatgcaatcacgatagatggcgtgactcgtttatcaacgcacgtcctccaaggcagtcacagctgtgcagagctgctctggagatgggtgccaagcccagctggagtcagcttccatccccaacgcacgtgagccgcccagtaaaaggtcaccatgaatgtcttgacat  
cctgatatcctgggcatagatgttgaccaagaattcctcatttgggaactcctctctatgtagcttctgtacagcaatccattgcactgtgaagcttctttatgctggtgctgactacagaaaggcaaatattgggatactccattacatgctgctgctcaacaatccagcacagaattgt  
aaacttactgctagaattggagcagatatcaatgccaaaatacagagcttctgcacatagatgtagctacgtctagcagatgggtgaaagaatatgtctcaacatgaagctaccccaagctctctttaccaactttgcccactctgtatccgaagctacataggaaaaccaagattgca  
ccttatcccaactccagctgccaagcttactgaagaatttctacagtatcgataa