

CST4-cystatin S Gene

Reference: PTXBC065714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cystatin-S; cystatin 4; cystatin-SA-III; salivary acidic protein 1; cystatin S		
For research use only.			

GENE DESCRIPTION	cystatin S
GENE ID	1472
ORIGIN SPECIES	Human
ACCESSIONS	BC065714
VECTOR	pDONR223
NCBI SYMBOL	CST4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgccctctgtgtacctgctactcctgatggctaccctggctgggctctggcctcagctccaaggaggagaataggataatccagggtggcatctatgatgcagacctcaatgatgagtggttacagcgtgcccttcacttcgcatcagcgagtacaacaaggccaccgaag
atgagtactacagacgcccctgcaggtgctgcgagccagggagcagaccttgggggggtgaattactcttcgacgtagaggtgggcgcaccatattaccaagtcgcccaacttggacacctgtgccctccatgaacagccagaactgcagaagaacagttgtgctcttcga
gatctacgaagtccctgggaggacagaatgtccctggtgaattccaggtgtcaagaagcctag