

VTCN1-V-set domain containing T cell activation inhibitor 1 Gene

Reference: PTXBC065717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B7-H4; B7H4; B7S1; B7X; B7h.5; PRO1291; VCTN1; V-set domain-containing T-cell activation inhibitor 1; B7 family member, H4; B7 homolog 4; B7 superfamily member 1; T cell costimulatory molecule B7x; immune costimulatory protein B7-H4; V-set domain containi		

For research use only.

GENE DESCRIPTION	V-set domain containing T cell activation inhibitor 1
GENE ID	79679
ORIGIN SPECIES	Human
ACCESSIONS	BC065717
VECTOR	pDONR223
NCBI SYMBOL	VTCN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcagagccgacagcagtggttctgatcaagtgatagttggcaatgcctcttggcgctgaaaaacgtgcaactcacagatgctggcacctacaaatgttatatcatcacttctaaaggcaagggaatgtaaccttgagtataaaactggagccttcagcatgccggaagtgaatgt
ggactataatgccagctcagagaccttgcggtgtgaggtctcccgatggttcccccagcccacagtggtctgggcatcccaagttgaccaggagccaaactctcggaagctccaataccagcttggagctgaactctgagaatgtgacctgaaggttggctgtcttacaatgttacg
atcaacaacatactctgtatgattgaaatgacattgccaaagcaacaggggatatacaagtgacagaatcggagatcaaaaggcgagtcacacagctgctaaactcaaaggcttctctgtgtctcttcttcttggccatcagctgggcacttctgcctcagcccttacctgatgc
taaaataa