

## C20orf70-chromosome 20 open reading frame 70 Gene

Reference: PTXBC065726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf70; PSP; SPLUNC2; bA49G10.1; BPI fold-containing family A member 2; parotid secretory protein; short palate, lung and nasal epithelium carcinoma-associated protein 2; BPI fold containing family A member 2		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	chromosome 20 open reading frame 70
<b>GENE ID</b>	140683
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC065726
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	C20orf70
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgcttcagcttggaaactgttctcctgtgcggcgtgctcactgggacctcagagctctcttgacaacttggcaatgacctaaagcaatgctgataagctggaacctgttcttcagaggacttgagacagttgacaatactcttaaaggcatccttgagaaactgaaggtcagccttagagtgctcagaaatcagtgcttggcaactggccaagcagaagcccagggaagctgagaaatgctgaacaatgtcatttetaagctgttccaactaacacggacattttgggtgaaaatcagcaactccctcatcctggatgtcaaagctgaaccgatgatgatggcaagccttaaactgagcttccctgtcaccggaatgtcactgtggccgggccatcattggccagattatcaacctgaaagcctccttggacctcctgaccgcagtcacaattgaaactgatccccagacacaccagcctgttgcgctcctgggagaatgcccagtgacccaaccagcatctccttctgctggacaacacagccaatcatcaacaagttcgtgaatagcgtgatcaacacgctgaaaagcactgtatcctccctgctgcagaaggagatattccactgatccgatcttccactccctggatgtgaatgtcattcagcaggtcgtc gataatcctcagcacaaccagcgtgcaaacctcatctga