

ZNF483-zinc finger protein 483 Gene

Reference: PTXBC065738

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZKSCAN16; ZSCAN48; zinc finger protein 483; zinc finger protein HIT-10; zinc finger protein with KRAB and SCAN domains 16		

For research use only.

GENE DESCRIPTION	zinc finger protein 483
GENE ID	158399
ORIGIN SPECIES	Human
ACCESSIONS	BC065738
VECTOR	pDONR223
NCBI SYMBOL	ZNF483
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagctgtagtgccttgaacaagatgacagccatccaccgaacctcaactctggcctgactgaacaaatgaggtcccaagagtggttactctggggaacaagaagctatttaagaggaaatgctgctgagcagagctttcagacagaggttaggtgtttgttactcaga
agtagctggaccaggaaagctctgagtcactctgggagctctgcaatcagtggtgagaccagacattcacacgaagaacagatttagagctctggtgtttgagcagttcctgaccattttgctggggagatcaggattgggtaaagtcacaacatcctgagagtagtgaggaa
ggtagccctaataagaatttgaccagatgctgaagaaaaagatccagtcctcagagattctactgtttcccaaggaggagactcaaaaggagataaaatggtcactgtttgtcccaatctgagtcctgtgaatctataacattgaagatgtagctgtgaactttcaagaggagagtgga
gaagctggagcctttcaaaaggagctatataaggagtgctactggaaacctcagggaacctagaatttctgactttccagtttcaaaattagagtgatttccagctaaagtgggtgaattgccatggctgctggaagaagctcaaaagctccgactaggagtggtattga