

## WDR90-WD repeat domain 90 Gene

Reference: PTXBC065836

|                |                                                                                                        |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C16orf15; C16orf16; C16orf17; C16orf18; C16orf19; WD repeat-containing protein 90; WD repeat domain 90 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | WD repeat domain 90                                                                               |
| <b>GENE ID</b>          | 197335                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC065836                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | WDR90                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctctgttggacagccccggcctggggctggaatggaatgcctgcgtccctttaccctgccaggggcccaacccccagcctagctacggcttccccaggtgaacgaggtggtcttcagccccggggagtcaccatgcgcccacatgcagtggatgggagtgctggg  
tggggccttggccagcatggagcttggatccagttccaggtgctgaaccagagctgctctgctggcatggagccccctgtgtggcgccctgagcagcagcggctagcggctggctacgggtgacggctccctgcgcctctcagcgtctcccgcacggccatggagctcaag  
atgcaacccccacccgggtggcgtgaccactgttgccttctccacccatgggtcagactgtcctctctggagacaagatgggctcgtggctgtgagccacccctgcacagggaacacttccgtgtgtgagtgaccaccaggggcgccecaatctctaccatctgtgcacgtgcaagaga  
ctcagggccacctgccaccctccctgccttgccttggatggggcgtcctgatgtacgtgggccccgggtgtttacaaggaggtgatcatctacaacctctgccagaagcaggtggtggagaagataccactgccctttttgccatgtccctgagcctgtccccgggacccac  
ctcctgctgttggcttctgagtgcatgctgaggtggtgagactgtccatggggactgccaagacttgcggccacgacaacgacgtgcacctgtgcaggtttacacctccgcagcgtgctcttcacggcccccgaacgagatccttgttgggaggtccccggcctctga