

CTDSP2-CTD (carboxy-terminal domain, RNA polymerase II, polypeptide A) small phosphatase 2 Gene

Reference: PTXBC065920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OS4; PSR2; SCP2; carboxy-terminal domain RNA polymerase II polypeptide A small phosphatase 2; CTD (carboxy-terminal domain, RNA polymerase II, polypeptide A) small phosphatase 2; NLI-interacting factor 2; conserved gene amplified in osteosarcoma; nuclear		

For research use only.

GENE DESCRIPTION	CTD (carboxy-terminal domain, RNA polymerase II, polypeptide A) small phosphatase 2
GENE ID	10106
ORIGIN SPECIES	Human
ACCESSIONS	BC065920
VECTOR	pDONR223
NCBI SYMBOL	CTDSP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagggctccatcatcaccaggcgcggagggaagacgcctgtgtctaccaagcaaggcctgtgtctccaagtcctctcctaagaagcctgtggacgtaacatctcaaggccctttctgtgttttgcgccagcattgtgccagtcaagttctccactgagctcgtgcgtataaggagggaagcaaacaccattgctaaagtcgagctgtctccagtgctcctcagttctaccagatcccaggacccctgtctccagaggtgacagagggaagatcaaggaaagatctgtgtgtcattgacctcgatgaacccctgtgcatagctctttaagccaatcaacaa
tgcctgacttcattagtgcttatagagattgaggggaccactcaccagggtatgtgtcctcaagggccttatgttgatgagttcctgagacgcattgggggaactctttgaatgtgttctctcactgccagcctggccaagtatgccaccctgtgacagacctgctggaccggtgtgggtgttcc
gggccgcctattccgtgagcttgcgtgttccaccagggtgctacgtcaaggacctcagccgctggggaggacgtgagaagaccctatcctggacaactgcctgttctacatattccaccccagagaatgcagtcctgtgcagtcctggtttgatgacatggcagacactgag
ttgctgaacctgatcccaatctttgaggagctgagcggagcagaggacgtctacaccagccttgggcagctgcggcccttag