

REEP5-receptor accessory protein 5 Gene

Reference: PTXBC065926

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf18; D5S346; DP1; TB2; YOP1; Yip2e; receptor expression-enhancing protein 5; deleted in polyposis 1; polyposis coli region hypothetical protein DP1; polyposis locus protein 1; receptor accessory protein 5		

For research use only.

GENE DESCRIPTION	receptor accessory protein 5
GENE ID	7905
ORIGIN SPECIES	Human
ACCESSIONS	BC065926
VECTOR	pDONR223
NCBI SYMBOL	REEP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagaggttcgaccgttctcgcacgagaagaactgcatgactgacctctggccaagctcgaggccaaaccggcgtgaacaggagcttcacgctcttgggtgcatcggactgggtggcctgtacgtggtgttcgggtatggagcctctctcctctgcaacctgataggattggct
accagcctacatctcaattaaagctatagagagtcaccaaaagaagatgataccagtggtgactactgggtagtgtatgggtgtgtcagcattgctgaattctctctgatactctctgtcatggttccctctactacatgctgaagtgtggcttctgttgggtgcatggccccgagcc
cttctaattgggctgaactgctcacaagcgcatcatcgtctcttctcctgaagcagcagtcagatggacagtgtgtgcaaggaccttaagacaaggccaaagagactgcagatgccatcactaaagaagcgaagaagctaccgtgaatttactgggtgaagaaagaagagcac
ctaa