

TMEM60-transmembrane protein 60 Gene

Reference: PTXBC065930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf35; DC32; transmembrane protein 60		

For research use only.

GENE DESCRIPTION	transmembrane protein 60
GENE ID	85025
ORIGIN SPECIES	Human
ACCESSIONS	BC065930
VECTOR	pDONR223
NCBI SYMBOL	TMEM60
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccttggctcagagagtactactcacctggctttcacactactcttcttgatcatgttggtgttgaacctggatgagaaagcaccttggaaactggctctcatattcattccagcttgatatttgatactatccttctgtcctgctgattgtgaaatggctggcggtgtaagtctggctttgacc
ctcgacatggatcacacaataattaaaaaaagcctgttacctcattgcaatgttacttaaatagcctctgcctcgcactctgtgtaaaactggaacagttactacatgaatctatcctatgtcttcattcctttatggccttgcctggctgggctttaacagaactcggatataatgtctttttgt
gagagactga