

NEK2-NIMA (never in mitosis gene a)-related kinase 2 Gene

Reference: PTXBC065932

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase Nek2; HsPK21; NEK2A; NLK1; PPP1R111; RP67; NIMA (never in mitosis gene a)-related kinase 2; nimA-like protein kinase 1; nimA-related protein kinase 2; protein phosphatase 1, regulatory subunit 111; NIMA related kinase 2		

For research use only.

GENE DESCRIPTION	NIMA (never in mitosis gene a)-related kinase 2
GENE ID	4751
ORIGIN SPECIES	Human
ACCESSIONS	BC065932
VECTOR	pDONR223
NCBI SYMBOL	NEK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccttccgggctgaggactatgaagtgtgtacaccattggcacaggctcctacggccgctgccagaagatccggaggagagtgatggcgaagatattgttggaaagaacttgactatggctccatgacagaagctgagaacagatgctgttctgaagtgaattgcttcgtgaac
tgaacatccaacatcgcttcttactatgatcggtatttaccggaccaatacaacactgtacattgtaattggaatatttgaaggaggatctggctagtgtaattacaagggaaccaaggaaaggcaatacttagatgaagagttgttctcagtgatgactcagttgactctggccct
gaaggatgccacagacgaagtgtgtgtgtatccgtatgcatcggtgatctgaaccagccaatgttttctggaagcgaagcaaaacgtcaagcttgagactttgggctagctagaatattaaccacgacacgagttttgcaaaaacattgttggcacacctattacatgtctctg
aacaatgaatgcctctacaatgagaatacagatctgtcattgggctgctgtatgagttatgtcattaatcttttag