

HIST1H2AJ-histone cluster 1, H2aj Gene

Reference: PTXBC066232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A/E; H2AFE; dJ160A22.4; histone cluster 1, H2aj; H2A histone family, member E; histone 1, H2aj; histone cluster 1 H2A family member j		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2aj
GENE ID	8331
ORIGIN SPECIES	Human
ACCESSIONS	BC066232
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2AJ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgctgtaagcaggaggcaagctcgcgccaaggccaagaccgctctctcggcgccggttcagttcccgtaggcccagtgcatcgccgtctccgcaaaaggcaactatcgaggcgggtcgggtctggagcgcgggtgtacctggcggcggtgctggagtacctg
accgccgagatcctggagctggctggcaacgcggcgccgcgacaacaagaagactcgcacatcccgctcaccctcagctggccatccgcaacgatgagagctcaacaagcttctgggcaaaagtcaccatgcacagggtggcgctcgtgcccaacatccaggccgtgctgcca
aagaaaactgagagccaccacaagactaagtaa