

HIST1H4I-histone cluster 1, H4i Gene

Reference: PTXBC066250

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/m; H4FM; H4M; histone H4; H4 histone family, member M; Histone 4 family, member M; histone 1, H4i; histone cluster 1, H4i; histone family member; histone cluster 1 H4 family member i		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4i
GENE ID	8294
ORIGIN SPECIES	Human
ACCESSIONS	BC066250
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4I
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgcggcaaggaggttaagggcctggggaagggggtgccaagcgccaccgcaaggtgctgcgcgacacatccaggggtatccaagccagccattcggcgccttgctcgcgcggcggtgaagcgcattctggcctcatatgaggagaccgcggagtgttgaaggtgttctggagaacgtgatccgggagcgcctgacctacacggagcagccaagcgcaagacggtcaccgcatggagctggtctacgcgtcaagcgccaggcgccgacacctatggttcggcgctaa