

IL21-interleukin 21 Gene

Reference: PTXBC066259

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CVID11; IL-21; Za11; interleukin-21; interleukin-21 isoform; interleukin 21		

For research use only.

GENE DESCRIPTION	interleukin 21
GENE ID	59067
ORIGIN SPECIES	Human
ACCESSIONS	BC066259
VECTOR	pDONR223
NCBI SYMBOL	IL21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagatccagtcctggcaacatggagaggattgcatctgtctgatgtcatctcttgggacactggtccacaatcaagctcccaaggctcaagatcggccatgattagaatgctcaactatagatattgtgatcagctgaaaaattatggaatgacttgtccctgaattctgccagc
tccagaagatgtagagacaaactgtgagtggtcagctttttctgctttcagaagggcccaactaaagtcagcaaatcaggaacaatgaaaggataatcaatgatcaatataaaagctgaaaggaaaccacctccacaatgcagggagaagacagaaacacagataacatgccctt
catgtattcttatgagaaaaaccaccaagaattcctagaagattcaaatcactctccaaaagatgattcatcagcatctgtcctctagaacacacggagtgagattcctga
```