

IL26-interleukin 26 Gene

Reference: PTXBC066270

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	AK155; IL-26; interleukin-26; interleukin 26		

For research use only.

GENE DESCRIPTION	interleukin 26
GENE ID	55801
ORIGIN SPECIES	Human
ACCESSIONS	BC066270
VECTOR	pDONR223
NCBI SYMBOL	IL26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgctggtgaatttcattttgaggtgtgggttctgttagtcaactctgtctcttgcattgcccaagcacaagcaatctctcaccaaaagttgtacccaaggggaacattgtccaagctgtgacgctctctatacaagcagcatggctcaaaagcaacgattcagaagaccgcataaaaa
atatagcattataaaaaaagaaacaaaaaacagttatgaaaaactgtcaattcaagaacagctctctcctctctcatggaagacgttttggtaactgcaattgcaaggctgcaagaaaaacgctttgtggaggactttcatagccttaggcagaattgagccactgtattcctgtgcttc
atcagctagagagatgaaatccattaccaggatgaaaagaatatttataggattggaacaagaagaaatctacaagccatcagtggaactggatattctcttctggattaaaaaattatggaaagcagtcagtaa
```