

IL9-interleukin 9 Gene

Reference: PTXBC066285

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP40; IL-9; P40; interleukin-9; T-cell growth factor p40; cytokine P40; homolog of mouse T cell and mast cell growth factor 40; p40 T-cell and mast cell growth factor; p40 cytokine; interleukin 9		

For research use only.

GENE DESCRIPTION	interleukin 9
GENE ID	3578
ORIGIN SPECIES	Human
ACCESSIONS	BC066285
VECTOR	pDONR223
NCBI SYMBOL	IL9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttctggccatggctcttacctctgcccctgctcctgtgtcctcgtggcaggccaggggtgtccaaccttggcggggatcctggacatcaactcctcatcaacaagatgcagggaagatccagctccaagtgccactgcagtgtctaatgtgaccagttgtctctgtttgggcattccctctgac
aactgcaccagaccatgcttcagtgagagactgtctcagatgaccaataccaccatgcaacaagataccactgatgttcagtcgggtgaaaaaatcagttgaagtactaaagaacaacaagtgccatattttctctgtgaacagccatgcaaccaaacccacggcaggcaacgcgtctgac
attctgaagagcttcttggaattttccagaaaaaagatgagaggatgagaggcaagatatga