

SPATA8-spermatogenesis associated 8 Gene

Reference: PTXBC066292

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SRG8; spermatogenesis-associated protein 8; spermatogenesis-related protein 8; testicular secretory protein Li 52; spermatogenesis associated 8		

For research use only.

GENE DESCRIPTION	spermatogenesis associated 8
GENE ID	145946
ORIGIN SPECIES	Human
ACCESSIONS	BC066292
VECTOR	pDONR223
NCBI SYMBOL	SPATA8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccccggctgggatgtcagggccaggacaactctgtctctaccagaaatgccctcttttcagaggtgccttgcctagaacctcgtcgcgacatttctcagaagccatgacatgtccctgcggctggaggccttcaagggtggccaggaggcctcaaggcccggtgtggcctgcaaggagaacagctgttctcacggaaggatccaaagggttcaaggcgaagggtgccatctgttcccactgattcagaagataaacaggagaagtgctctgtttcaccatattgctggtctga