

LYZL2-lysozyme-like 2 Gene

Reference: PTXBC066294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LYZD2; lysozyme-like protein 2; lysozyme D2; lysozyme-2; lysozyme like 2		

For research use only.

GENE DESCRIPTION	lysozyme-like 2
GENE ID	119180
ORIGIN SPECIES	Human
ACCESSIONS	BC066294
VECTOR	pDONR223
NCBI SYMBOL	LYZL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggacgtccctgagctgctgacccgactaagtggagcagtggtttctccgagactcaactgagaagtcagcctctgcggcaggcaccagggaatctgctttcagttctgtctccggcaggtttgaggatgaaggctgcgggcatttgaccctcattggctgcctgtcaca
ggcgcgagtcacaaatctacactcgttgcaactggcaaaatattctcagggtgctggaacttactgggcttcagccttggaactggatctgcatggcgtattatgagagcggctacacaccacagcccagacggctcgtgatgacggcagcatcactacggcatctcc
agatcaacagcttcgctggtgagacgcggaagctgaaggagaacaaccactgccacgtcgcctgctcagccttggtcactgatgacctcacagatgcgattatctgtgccaagaaattgttaaagagacacaagggaatgaattatggcaaggctggaagaacactgtgagggga
gagacctgtccgactggaaaaagactgtgaggtttcctaa