

ATP5J-ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit F6 Gene

Reference: PTXBC066310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP5; ATP5A; ATPM; CF6; ATP synthase-coupling factor 6, mitochondrial; ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit F6; ATPase subunit F6; coupling factor 6; mitochondrial ATP synthase, coupling factor 6; mitochondrial ATP synthase, su		

For research use only.

GENE DESCRIPTION	ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit F6
GENE ID	522
ORIGIN SPECIES	Human
ACCESSIONS	BC066310
VECTOR	pDONR223
NCBI SYMBOL	ATP5J
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattcttcagaggctcttcaggtctctctgtcattcgggtcagccgtctcagtcatttgcggagggaacattgggtgttacagcagtggtcatttaataaggaaacttgatcctatacagaaactcttgggacaagattagagaatacaaaactaagcgacagacatctggaggacctgttgatc
tagttcagagtatcagcaagagctggagagggagcttttaagctcaagcaaatgttggtaatgcagacatgaatacatttcacacctcaaatgtgaagatcccaaatgtgaagtcacgaaaaacccagcctga