

XAGE2-X antigen family, member 2 Gene

Reference: PTXBC066311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT12.2; GAGED3; XAGE-2; XAGE2B; X antigen family member 2; G antigen family D member 3; X antigen family, member 2B; cancer/testis antigen 12.2; cancer/testis antigen family 12, member 2		

For research use only.

GENE DESCRIPTION	X antigen family, member 2
GENE ID	9502
ORIGIN SPECIES	Human
ACCESSIONS	BC066311
VECTOR	pDONR223
NCBI SYMBOL	XAGE2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgaggagatcaacatagggcctaggccaagaaggtttacagcctctgagctgattgggctatgctgaaccactgatgaagagcctaagaagagaaccacccactaaaagtcggaatcctacacgtatcagaagagagaagatgacagggtgcagctga
gattcaagtgcctgacctggaagccgatctccaggagctatgtcagacaaagactgggagtgatgtgaaggtgtgactgatgtcaaggggaagattctacaaaagcagagcactttaaatgccagaagcaggtgaagggaatcacaggtttaa