

RPL18A-ribosomal protein L18a Gene

Reference: PTXBC066319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L18A; 60S ribosomal protein L18a; ribosomal protein L18a-like protein		
For research use only.			

GENE DESCRIPTION	ribosomal protein L18a
GENE ID	6142
ORIGIN SPECIES	Human
ACCESSIONS	BC066319
VECTOR	pDONR223
NCBI SYMBOL	RPL18A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggcctcgggcacgctacgagagtaaggtagtggtcgtgcctgccacccccaaatgccacacgcgcctctaccgcatgcgaatcttgcgcctaatcatgtcgtcgccaagtcccgcttctggtactttgtatctcagttaaagaagatgaagaagtcttcaggggagatt
gtctactgtggcaggtgtttgagaagtcctcctcgggtgaagaacttcgggatctgctgcgtatgactcccgagcggcaccacaacatgtaccgggaataccgggacctgaccacagcaggcgtgtcaccagtgctaccgagacatgggtgccggcaccgcgccga
gcccactccattcagatcatgaagggtggaggatcgcgcgcaagtgccgcggcggcgtgtcaagcagttccagactccaagatcaagttcccgctgccccaccgggtcctgcgcgtcagcacaagccacgttcaccaccaagaggccaacaccttcttag