

RPL17-ribosomal protein L17 Gene

Reference: PTXBC066324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L17; PD-1; RPL23; 60S ribosomal protein L17; 60S ribosomal protein L23; gene encoding putative NFkB activating protein; ribosomal protein L17		

For research use only.

GENE DESCRIPTION	ribosomal protein L17
GENE ID	6139
ORIGIN SPECIES	Human
ACCESSIONS	BC066324
VECTOR	pDONR223
NCBI SYMBOL	RPL17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttcgtctattcacttgacccggagaacccacgaaatcatgcaaatcaagaggttcacatctcgtgttcactttaagaacactcgtgaaactgctcagggcatcaagggtatgcataacgaaagccacgaagtatctgaaagatgtcactttacagaacagtggtaccattccgacg
ttacaatggtggagttggcaggtgtgcgcagggccaagcaatgggctggacacaaggtcgggtggcccaaaagagtgctgaattttgctgcacatgcttaaaacgcagagagtaatgctgaactaaagggttagatgtagattctctggtcattgagcatatccaagtgaacaaagcacc
taagatgcgcgcggacacacagagctcgtggtcggattaaccatacatgagctctccctgccacattgagatgatccttacggaaaaggaacagattgtcctaaccagaagaggaggtgcccagaagaaaaagatatccagaagaactgaagaacaaaaacttatggcagc
ggagtaa