

PSMD14-proteasome (prosome, macropain) 26S subunit, non-ATPase, 14 Gene

Reference: PTXBC066336

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAD1; POH1; RPN11; 26S proteasome non-ATPase regulatory subunit 14; 26S proteasome regulatory subunit rpn11; 26S proteasome-associated PAD1 homolog 1; proteasome (prosome, macropain) 26S subunit, non-ATPase, 14; testis tissue sperm-binding protein Li 69n;		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) 26S subunit, non-ATPase, 14
GENE ID	10213
ORIGIN SPECIES	Human
ACCESSIONS	BC066336
VECTOR	pDONR223
NCBI SYMBOL	PSMD14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagactcttagactggaggaggtatgctggagggccagggccacctacagatgctcctgcagtggaacacagacagaacaagtctatatctctccctggcactgttaaaatgttaaaacatggccgtgctggaggtccaatggaagtattgggttgatgcttgagaattgt
tgaatgattataccgtcagagtgattgatgtttgctatgccacagtcaggacaggtgtcagtggtggagggcagttgatccagtggtccaagctaaatgttgatattgtgaagcagacaggaagccggagatggtgttggttggtatcacagtcacctggccttggttggtctctggt
gtggatatacaactcagcagagctttgaagcctgtcggagagagctgtggcagtggtgtggatccatcagagtgtaaaaggaaaggtgttattgatccctcagattgatcaatgctaataatgatgtcttaggacatgaaccaagacaacaactcgaatctgggtcacttaacaag
ccatctatccaggcattatcatggactaaacagacattactcattactattaactatcgaaaaatgaactggaacagagaagatgttctaattgtcataagaagagttggatggaaggttgacacttcagcactacagtgacattgtaaacacaatgaatcagtggttaaggagatgtt
ggaattagccaagaattacaataaggctgtagaagaagaagataagatgacacctgaacagctggcaataaagaatgttggcagcaggaccccaaacgtcatttggaggaaacatgtggtgacttatgacctcaaatattgtccagtggttagcagctatgttgatactgtctatttaaat
aa