



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C20orf71-chromosome 20 open reading frame 71 Gene

Reference: PTXBC066354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf71; SPLUNC3; BPI fold-containing family A member 3; short long palate, lung and nasal epithelium carcinoma associated 3; BPI fold containing family A member 3		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 71
GENE ID	128861
ORIGIN SPECIES	Human
ACCESSIONS	BC066354
VECTOR	pDONR223
NCBI SYMBOL	C20orf71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatagtcttcaccttggaggctctcattctctcggggtgctggcccttggccacacaaagcagcttggctgctggcccaagcccaagagacacaaatccaccctggcaagaattatgtctcaggggctctataagcacaaacgacgaagccgaatfcagaacatcc
actttggggacagactgaattgctcagcacaaatgggccccagggtctgggggctgctaatcagcggcaggaacaccacagcagcagaagagagcaggtctctcagataacacatgttaaagatgtgtgcacatatgagcatgtttggagttctggctggagaaagacagattggc
cggggagatctgtgataggcaaaftcgatfcagagcccgactgtctcattggtctcattctcactgaggtctcaccaccaaagatgaatcagttttctacacctcaagaagaatctgcaaaaagttctccacacatggtgaagctcaggtatgcttcctgacgttccccaactgta
cagctgattgtgaaaacttgaaaagctcagacagagagctctcattgaaccacacccacatccacacagctcgtcagtcgacgaagctcgcacagctccccaactgta