



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SIN3A-SIN3 homolog A, transcription regulator (yeast) Gene

Reference: PTXBC066364

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

transcriptional regulator, SIN3A; transcriptional corepressor Sin3a; transcriptional corepressor Sin3A; histone deacetylase complex subunit Sin3a; paired amphipathic helix protein Sin3a; WITKOS; SIN3 homolog A, transcription regulator; SIN3 transcription

For research use only.

GENE DESCRIPTION	SIN3 homolog A, transcription regulator (yeast)
GENE ID	25942
ORIGIN SPECIES	Human
ACCESSIONS	BC066364
VECTOR	pDONR223
NCBI SYMBOL	SIN3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcggcggttgatgaccaggaggtacccgggtgatgcagccagcagcgtcggatccctggcagcacagagggctttctaccagaccggggtcgtccccctgccctctgtgtatgaagcagtgctctgagaccatgcagtcagctacgggaattcagtactctgtaacacca
gctacaggttgatcagctccagcagcagctccggcgatctcggcgccgctatagcagcagtcattacagcattacacccaacagcggcgtagcccccagggcgagcagctgtccagagctatgcagctacccacagctgcagcagctcagggcagcagacagcattcagag
ctcgaaggttgatcagctcgtctttcttggctgcaaaaatttctatgcaggaataatctgtgtcagctgtgtgtgtgtgtacatacttcaactcagagatactag

More info: Please contact us for any question: contact@proteogenix.fr