

NCRNA00052-non-protein coding RNA 52 Gene

Reference: PTXBC066551

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00052; TMEM83; long intergenic non-protein coding RNA 52		
For research use only.			

GENE DESCRIPTION	non-protein coding RNA 52
GENE ID	145978
ORIGIN SPECIES	Human
ACCESSIONS	BC066551
VECTOR	pDONR223
NCBI SYMBOL	NCRNA00052
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattgtgatgcttgcacatcattctgcaatccagaagatttttgcataatttttgcattacagaaaaatctcagtcctccctctctctctcaatctgtgtgtctcttttactccatatctctgtgtgtgtcttttactccatatctctgtgtgtgtgtttatgtctctctctctc
cttcccatgtttctctctcacacacacacactcattcacagcttcaaaagacacgtctgtcttaccttcactttttgttttaa