



Reference: PTXBC066594

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC100126583
GENE ID	100126583
ORIGIN SPECIES	Human
ACCESSIONS	BC066594
VECTOR	pDONR223
NCBI SYMBOL	LOC100126583
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]