

C9orf37-chromosome 9 open reading frame 37 Gene

Reference: PTXBC066646

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C9orf37; AD038; chromosome 9 open reading frame 37		
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 37
GENE ID	85026
ORIGIN SPECIES	Human
ACCESSIONS	BC066646
VECTOR	pDONR223
NCBI SYMBOL	C9orf37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttactatgttcccaggtgatctggaactcctgacctcaagcaatccacctgcctcagctcccaagtaactgggattacaggaggagccaccgtgccggcctggctctgtgcattttatagataaggtaactgacaaacctagccacagtcaccccttgccttaaggag
aactggaattgaaccagagccttctcaccaccatctcgttgtttctggaagcaccctccaggcaggccagccagcatcacgggctccctggggctggaagctctgcaggatgccctttgagaaatgctgttcacagaaacctgaccttcaggccttgagacgtgggcagg
gagaagcagcgtccctcagagccagcctggcagtggtgctagcaggggccaaggccaggaggagcagggtctcctgtcggaggacgtgggcaagccctccacgcgcagcgggtttctcagcagggagggtccacaccacaccgttgggaacctgggtgcctaa