

DSCR9-Down syndrome critical region gene 9 (non-protein coding) Gene

Reference: PTXBC066653

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00038; Down syndrome critical region gene 9 (non-protein coding)		
For research use only.			

GENE DESCRIPTION	Down syndrome critical region gene 9 (non-protein coding)
GENE ID	257203
ORIGIN SPECIES	Human
ACCESSIONS	BC066653
VECTOR	pDONR223
NCBI SYMBOL	DSCR9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcaggatttccccgtgaactcccgccgacgcaggtccgcgcgagggccggcgccgagcgggcattccctccctatcacagctccaggggtgtgccccccggctgtggagccagaccccgccgtccagcagcgtatagacgagcccatgtctgcgatgtaactgc
accgcctggggatccacctcccgccagggtgagggcgccgtctccagcgtatgctaggcagcgcgcccgccctctgtgacgtgaccacgcctacagccgaaggatgggaaatgcggccctgttctgaacaggaagcggagcaggaaggaactgagaacaccag
catctgtccgcccgtgaatcagcgttagcatgtgaggcaacctcagccccagttttgaagaagggtgattcagattcagacacaggaattatggtga