

MRFAP1L1-Morf4 family associated protein 1-like 1 Gene

Reference: PTXBC066897

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PP784; MORF4 family-associated protein 1-like 1; Morf4 family associated protein 1 like 1		

For research use only.

GENE DESCRIPTION	Morf4 family associated protein 1-like 1
GENE ID	114932
ORIGIN SPECIES	Human
ACCESSIONS	BC066897
VECTOR	pDONR223
NCBI SYMBOL	MRFAP1L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggccctggacatagacgaggtggaagcgcctgagggaagtggaggtgctggagcccaggaggatttcgagcagttcctgctcccgtcatccacgagatgcgcgaggatcgcgtcttatacgcgagcacggcggtgtacctgctgaccaggagcaagctgtggga
gatggacaatatgcttatccagatcaaaacgcaggtggaggcctggaggagagcgcctcaatcacgtgcagcacccgagtggcgaagccgacgagagagtgctcggagttgtgcgagaagcgtgaggagaaagccaaggagattgcgaagatggcagagatgctgtcgagctc
gtctggcgaatagagagaagcgagctcttga