

RAB34-RAB34, member RAS oncogene family Gene

Reference: PTXBC066904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB34, member RAS oncogene family; RAB39; RAH; ras-related protein Rab-34; ras-related protein Rab-39; ras-related protein Rah		

For research use only.

GENE DESCRIPTION	RAB34, member RAS oncogene family
GENE ID	83871
ORIGIN SPECIES	Human
ACCESSIONS	BC066904
VECTOR	pDONR223
NCBI SYMBOL	RAB34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacattctggcaccctgtcggaggatcgcgtctgctggcggagctgccccagtgctgaggaaggagccgcttgcacgggcacaaagacttccacccccgcgtcacctgcgcctgccaggagcaccggacaggcaccgtgggatttaagatctcaaggtcattgtgtgggg
gacctgtcgggtgggaagacttgcctcattaataggttctgcaaaagacaccttfgataagaattacaaggccaccattggagtggaacttcgagatggaaacgatttgagtgctggcattcccttcagtttgagctttgggataccgctgggcaggagaggttcaaatgcatcaaccta
ctatagaggagctcaagccatcatcattgtcttcaacctgaatgatgtggcatctctggaacataccaagcagtgctggccgatgccctgaaggagaatgaccttccagtggtcttcttctttaggttccaagaagatctgagtagccctgctcagtagcgtgatggagaagacg
ccctccaggtggccaggagatgaaggctgagtagtggcagtcctctcactggtgagaatgtccgagaattcttctcgtgtggcagcactgaccttggagccaatgtgctggctgagctggagaatcgggggtcgcgcattggggtgttgcgcacacagtgatgac
aacaacctctacctaactgccagcaagaagaagcccacatgttggccatga