



Reference: PTXBC066909

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIAR; ankyrin repeat domain-containing protein 54; lyn-interacting ankyrin repeat protein; ankyrin repeat domain 54		

For research use only.

GENE DESCRIPTION	ankyrin repeat domain 54
GENE ID	129138
ORIGIN SPECIES	Human
ACCESSIONS	BC066909
VECTOR	pDONR223
NCBI SYMBOL	ANKRD54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1