

PPP1R2P3-protein phosphatase 1, regulatory (inhibitor) subunit 2 pseudogene 3 Gene

Reference: PTXBC066922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	protein phosphatase 1, regulatory (inhibitor) subunit 2 pseudogene 3
GENE ID	153743
ORIGIN SPECIES	Human
ACCESSIONS	BC066922
VECTOR	pDONR223
NCBI SYMBOL	PPP1R2P3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctcgacggcctccaccggcccatcaaggggatcttgaagaacaagacctctacgactctctatggtggcgctcgccgaacagccccgcaggagtgtcgacgaggagctgagcaaaaaatccagaagtgggatgaaattaacatcttggcgacatcatccagcag
acaaaggctatggttaatgaaaatagatgaaccaagccctctaccatagatgatgggtgatgaagatgcgtgttagggacaccgagaccactgaagccatggcgccaggcatcctagccaaagaattagctgctgctgaaggcttgagccaaagtaccggattcagggaacaag
aaagcagtgagaggagatagtgacctctcactgaagaacgagaaaaaaacgacaatttgaatgagaagggaagcttcactacaatgaaggactcaatatcaaaactagccagacaattaatttcaaaagacctacatgatgatgaagatgaagaaatgttagagactgcagatgg
agaaagcatgaatacggagaatcaaatcaagatctactccaagtgaaccaacagcaaaacaattacgaagttcatag