

SETD3-SET domain containing 3 Gene

Reference: PTXBC066931

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone-lysine N-methyltransferase setd3; C14orf154; SET domain-containing protein 3; SET domain containing 3		

For research use only.

GENE DESCRIPTION	SET domain containing 3
GENE ID	84193
ORIGIN SPECIES	Human
ACCESSIONS	BC066931
VECTOR	pDONR223
NCBI SYMBOL	SETD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtaagaagagtcgagtaaaactcagaatctggcactggctacagcaactgtgtcaccaggaatcttgaaactgaccagtgagctgtgcagaaatgcagcagtcggcgccctggccaggaagagtgaggagagatgtgcagatccggactctggttgagaaat
acggaaaaagcaaaaaggctgtcgttacttttgatggaaaaagagaagattacttctgatctaataaaggcccttgaaatgggcttctgtcaggggtttgaaatggttaactcaaaagagaggccttgggttgagagcaacaagagatatcaaggcagaagaattgtttatg
gggtccacgaaatgtgaatgactgtgaatctgctaaaaatcagtggtggggcccttatattctcaagaccgaatcctcaagccatgggaacatcgcactggcctttcattgtgtgtgagcagccagccctaacctctctggcagccctatatcaaacctcccagtgaaatgaca
ctctctctactttgaagaagatgaagttcggatcttcagtcacacaaagctatacatgatgtcttcagccagtatataaaacacagctcgcacagtagcctacttataaagtcacagacccatctcatgccaacaactacccttgaggattctttcactacgaggactacaggtgggc
agtctctctgttatgacgaggcaaaaccaaattcccacagaggatggttcccgctgacccctggctctgattcctttatggatatgttaaccacacaaacggcctggaacacctcttcccctgtgtgtgtga