

STAG3L3-stromal antigen 3-like 3 Gene

Reference: PTXBC066937

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STAG3L1; STAG3L2; OTTHUMP00000211421; OTTHUMP00000211422; OTTHUMP00000211424; STAG3-like protein 3; stromal antigen 3-like protein 3; stromal antigen 3-like 3		

For research use only.

GENE DESCRIPTION	stromal antigen 3-like 3
GENE ID	442578
ORIGIN SPECIES	Human
ACCESSIONS	BC066937
VECTOR	pDONR223
NCBI SYMBOL	STAG3L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatctttcaatgctgcgaaagctgccaaaagtacatgcagggatgctcttctgagatccgtgctatctgcattgaggaaattgggtgtggatgcaagctacagcacgtcttctcaccgacagctatttaaaatatattgggtgactctgcatgataagcaccgagaagtcgcgtga
agtgcgtgaaggctctgaaagggtgtacgttaaccgggacctgaccgcacgtggagctcttactggccgttcaaggacgtggtttccatgatcatggacagagagtacagtgtggcagtgaggccgtcagattactgatacttatccttaagaacatggaaggggtgctgat
ggacgtggactgtgagagcgtctacccattgtgtag